

ANALYTICAL PERFORMANCE ASPECTS OF F-POINT FIBRINOGEN IVD ASSAY

ACCURACY, PRECISION & LIMITS OF DETECTION



TASKS

- I. Within-run Accuracy and Precision
- II. Between-run Accuracy and Precision
- III. Limit of Detection (LoD) and Lower Limit of Quantification (LLOQ)

GUIDANCE

"Guideline on bioanalytical method validation". EMEA/CHMP/EWP/192217/2009 Rev. 1 Corr. 2. – European Medicines Agency, 2011.

ANALYTICAL SAMPLES:

Analytical standards are made from the set of four artificial fluids using the following mixing formula:

$$(M:B:P)+S$$

where M, B, P are the volumes of corresponding fluids and S is amount of final fibrinogen concentration spike (g/L)

M	• Matrix: The Model 046, Blood Mimicking Fluid (CIRS Inc.);
B	• Buffer: Imidazole (HYPHEN Biomed)
P	• 222101 Plasma Calibrator (HYPHEN Biomed)
S	• Spike: Fibrinogen, Plasminogen-Depleted, Human Plasma (95% clottable, MERCK Inc.)



ACCURACY AND PRECISION

QUALITY CONTROL PROCEDURE

4 Quality Control (QC) Samples -> Number of runs
 -> F-Point Response -> Back-calc. concentrations
 -> Accuracy as mean difference with nominal values & Precision as CV

QC levels	Prep Protocol {M:B} (v/v) + S (g/L)	
QC **	0.81	{0.2:0.8} + 1.0
QC Low	1.59	{0.2:0.8} + 2.0
QC Med	3.17	{0.2:0.8} + 4.0
QC High	4.72	{0.2:0.8} + 6.0

BETWEEN-RUN

Quality Control	conc X		back-calc conc X	Accuracy, %	mean % diff	CV%	acceptance criteria, +/- %
QC **	0.81	Day1	0.9	119.3	19.3	9.8	20.0
		Day2	0.9				
		Day3	1.1				
QC Low	1.59	Day1	2.1	113.2	13.2	12.0	15.0
		Day2	1.7				
		Day3	1.6				
QC Med	3.17	Day1	2.9	95.7	-4.3	3.1	15.0
		Day2	3.1				
		Day3	3.1				
QC High	4.72	Day1	4.6	98.9	-1.1	2.0	15.0
		Day2	4.8				
		Day3	4.6				

BETWEEN-RUN
QUALITY CONTROL:
Passed

WITHIN-RUN

Quality Control	conc X		back-calc conc X	Accuracy, %	mean % diff	CV%	acceptance criteria, +/- %
QC **	0.81	Replicate1	0.9	116.0	16.0	15.9	20.0
		Replicate2	0.8				
		Replicate3	0.8				
		Replicate4	1				
		Replicate5	1.2				
QC Low	1.59	Replicate1	2.1	110.7	10.7	10.5	15.0
		Replicate2	1.6				
		Replicate3	1.6				
		Replicate4	1.7				
		Replicate5	1.8				
QC Med	3.17	Replicate1	2.9	90.2	-9.8	4.7	15.0
		Replicate2	3				
		Replicate3	2.9				
		Replicate4	2.6				
		Replicate5	2.9				
QC High	4.72	Replicate1	4.6	106.4	6.4	4.8	15.0
		Replicate2	5				
		Replicate3	5.2				
		Replicate4	5.3				
		Replicate5	5				

WITHIN-RUN
QUALITY CONTROL:
Passed



LIMIT OF DETECTION & LOWER LIMIT OF QUANTIFICATION

LOD

6 blank Samples -> F-Point -> Mean and noise evaluation

-> $LoD = 3.3 * St.Dev + Mean$

	CS sample	fibrinogen, g/L	back-calc. value, g/L
1	Blank	0	0.00
2	Blank	0	0.15
3	Blank	0	0.10
4	Blank	0	0.00
5	Blank	0	0.15
6	Blank	0	0.00
mean :			0.066
stdev :			0.069
LoD :			0.293

LOD:
0.29 g/L

LLOQ

3 Candidates: LoD, 2*LoD, lowest QC -> F-Point

-> Assessment of Accuracy & Precision

-> LLoQ = min conc where acceptance criteria met

Candidate Sample	conc x	back-calc. value, g/L		mean% diff	CV%	Acceptance criteria%, +/-
LoD	0.29	within-run	5 runs	31.0	29.4	20
		between-run	3 days	26.1	11.7	20
2*LoD	0.59	within-run	5 runs	15.3	21.1	20
		between-run	3 days	19.1	10.5	20
Lowest QC	0.81	within-run	5 runs	16.0	15.9	20
		between-run	3 days	19.3	9.80	20

LLOQ:
0.81 g/L